

# PRACTICAL NO: 1

Date : 09/01/2025

TITLE: To simulate Electroencephalogram Signal.

AIM / OBJECTIVE: To understand normal size EEG pattern and various abnormalities associated with EEG.

APPARATUS / TOOLS / EQUIPMENT / RESOURCES USED:

Biomedical and Signal processing lab virtual lab.

PC

Reference book.

Notes.

CONCEPT / THEORY OF EXPERIMENT:

The EEG is a unique and valuable measure of the brain's electrical functions. It is a graphic display of a difference in voltages from two sides of brain function recorded over time. EEG involves the study of recording these electrical signals that are generated by the brain.

PROCEDURE:

1. Select the type of patient from the waveform Selector menu.
2. Set the desired sampling rate and the number of samples to be displayed on screen.
3. Run the simulation.
4. Observe the various waveforms at different lead as shown in the fig.
5. Go to simulate tab for performing the simulation.

Parameter	F7	F3	FP1	FP2	F4	F8	T3	C3	C4	T4	T5	P3	P4	T6	CZ	O1	Fz	O2	Pz
Healthy Patient Baseline Activity	73.6	174	277	30.4	38	-7.71	138	232	19.5	54.9	54.7	163	130	235	45.1	104	85	64	
Left hand Backward	182.6	44.6	338	-667	194	-234	-4.2	59.4	-12.9	-286	22.5	123	123	73.9	19.6	142	187	79.5	
Left hand Forward	167	148	172	-267	157	-393	-37.2	81.9	-16.6	-98.3	197	119	119	7.5	87.2	119	192	89.2	
Left hand Imaginary Backward	151	446	264	-116	210	-221	-14.5	74.9	10.6	117	58.7	90.6	90.6	77	136	134	261	71.8	
Left hand Imaginary Forward	167	148	172	-267	157	-393	-37.2	81.9	-16.9	-99.2	203	119	119	7.32	87.6	119	192	98.6	
Left hand Random Movement	189	56.7	145	-143	88.5	192	224	110	-11.9	57.2	26.6	-50.5	-50.5	60.3	113	79.7	40.1	12.9	
Left Leg	144	228	172	-120	200	-199	36.5	21.8	38	-165	41.2	110	110	103	-79.7	100	198	55.4	
Right hand Backward	130	238	-26.9	-59.6	116	-249	-24.4	90.5	24.9	-59.7	-53.5	111	111	38.2	256	82.7	168	107	
Right hand Forward	132	269	91.7	-317	61.6	-45.7	17.9	59.4	54.9	46.1	29.9	96.7	59.3	266	57.9	151	104	106	
Right hand Imaginary Backward	207	199	144	-255	141	-405	26.7	75.6	-12.9	-138	63.3	121	121	42	-74	111	237	100	
Right hand Imaginary Forward	132	212	12.6	-331	-31.9	-428	-22	65.6	4.3	-354	-401	89.5	89.5	59.5	78	99.1	187	83.9	
Right hand Random Movement	158	148	172	-273	157	-393	-37.2	81.9	-16.6	-99.2	203	119	119	10.6	87.6	119	192	96.7	
Right Leg	84.4	131	124	-50.6	208	-121	34.2	53.4	-3.13	-69.7	42.8	83.1	83.1	18.2	38.3	95.8	83.8	113	

Table of election Readings.

## OBSERVATIONS

### ACTION POTENTIAL.

The information transmitted by a nerve is called an action potential APs are called by an exchange of ions across the neurons membrane and an AP is an temporary change the membrane potential that is transmitted along the axon. It is usually initiated in cell body and normally travels in one direction.

### Steps in AP Generation in Human Brain.

1. Nerve cells dendrites receive stimulate and Na<sup>+</sup> channels channels opens.
2. Action threshold is reached, causing additional Na<sup>+</sup> channels to open and depolarization occurs driving the membrane potential upto approximately +30mV.
3. Na<sup>+</sup> channels close and K<sup>+</sup> channels open, leading to repolarization back towards the rest potential.
4. Repolarization overshoots to about 90mV caused by hyperpolarization which prevents new stimuli and ensures signal transmission in one direction.

## CALCULATIONS :

### BRAIN RHYTHMS.

There are 5 major brain waves distinguished by their freq range.

Delta : below 5Hz

Theta : 4Hz to 8Hz

Beta : Above 13Hz

Gamma : 30 to 40 Hz

Alpha : 8 to 13Hz

## RESULTS :

- Observed the normal EEG pattern
- Understand other abnormalities for various associated with EEG.
- Observed the EEG waveform for various different body movements

## CONCLUSION :

Simulated EEG signal and observed the different abnormalities associated with EEG.

## PRACTICAL NO: 2

Date : 16/01/2025

TITLE: To simulate Electrocardiogram waveform.

AIM / OBJECTIVE: To understand normal ECG waveform of also various abnormalities associated with ECG.

APPARATUS / TOOLS / EQUIPMENT / RESOURCES USED:

Biomedical and signal processing lab - Virtual lab.

Pc







Reference book.


CONCEPT / THEORY OF EXPERIMENT:

Electrocardiogram signal is a graphical representation of the electrical activity occurring in the heart over time. The heart's electrical activity is generated by specialized cells within the heart's conduction system. This electrical activity coordinates the heart muscles contraction and relaxation.



PROCEDURE :

1. Open the simulation page on Virtual lab.
2. Select the appropriate ECG waveform.
3. Use the ON OFF toggle to stop the waveform and label the different.
4. Simulate the ECG waveform for different heart rate and save the output for each.
5. Compare and analyze the results for normal and abnormal waveform.





(C)



## OBSERVATIONS

### ECG Generation:

The human heart's electrical system comprises of 3 main parts. The Sinoatrial (SA) node, a trionodal (AV) node and the His-Purkinje system. Each heartbeat involves diastole (relaxation) and systole (contraction) phases of the heart beat, acting as the natural pacemaker, signaling the atria to contract and fill the ventricles with blood. The signal then travels to the AV node, which briefly delays, it allowing the ventricles to fill completely. Subsequently, the signal spreads through the His-Purkinje system, causing coordinated ventricular contraction. This electrical activity can be measured by electrodes placed on the body surface, recording an ECG.

### QRS complex

Description - Second major deflection seen on the ECG. Significance

Significance - Represents ventricular depolarization.

Duration - Naturally 60 to 100 milliseconds.

Shape - Can vary in morphology depending on lead and orientation of the heart.

### P wave

Description - First deflection seen on ECG.

Significance - Represents ventricular depolarization.

Duration - 60 to 100 milliseconds.

Shape - Smooth and rounded.


### T wave

Description - Follows the QRS complex.

Significance - Represents time taken for electrical impulse to travel from atria to ventricles.


Duration - About 120 to 100 milliseconds.

## CALCULATIONS



*Structure of Heart*


## RESULTS :



*Circulatory System.*

## CONCLUSION :

Observed Various abnormalities in ECG waveforms at different heart rates.



Electromyogram  
Signal.

## PRACTICAL NO: 3

Date : 30/01/2025

TITLE : To Simulate Electromyogram signal.

AIM / OBJECTIVE: To simulate EMG waveform and analyze the same.

APPARATUS / TOOLS / EQUIPMENT / RESOURCES USED:

Biomedical and Signal Processing Virtual Labs.  
Internet.

Reference Book.

Pc.

CONCEPT / THEORY OF EXPERIMENT:

Electromyography is a technique for evaluating and recording the electrical activity produced by skeletal muscle. EMG is performed using an instrument called an electromyograph to produce a record called an electromyogram.

PROCEDURE :

1. Open the Simulation Tab on the virtual lab.
2. Generate Gaussian random signals by changing the frequency.
3. Observe the deviation gaussian function by various contols.
4. Truncate the signal by changing the alpha.
5. Observe the power spectrum of EMG.

## OBSERVATIONS

### \* Electromyography.

The evolution and records the electrical activity of the skeletal muscles using an EMG. This technique detects muscle cell electrical potential when activated. Two common types are intramuscular. EMG using needle electrodes inserted into muscle tissues, and surface EMG (SEMG), recording muscle activity from the skin surface. SEMG is non invasive and widely used for diagnosing neuromuscular disorders, rehabilitation and device control application.

### \* EMG Measurement.

The human muscular system is composed of muscular fibre organised into motor units, each controlled units, each controlled by the single motor neuron. Muscle contraction, essential for movement, involves types of Isometric (muscle length remains constant); Concentric (muscle remains shortens) and Eccentric (muscle lengthens). Motor units are recruited based on the size principle. Smaller units are recruited based on the size principle and larger units as needed for the increased force. Motor units action potentials (MUAPs) propagate along muscle fibres, generating electrical activity detected, as fibre arrangement, electrode proximity and tissue properties. Surface EMG record microvolt level signals. Motor unit fire at varying rates, influencing muscle fibres type and function.

## CALCULATIONS :

Reading 2

Formula for EMG:

$$\text{Velocity } v_1 = \frac{l_1 - l_2}{t_1 - t_2}$$

$$v_1 = \frac{l_1 - l_2}{t_1 - t_2}$$

From Virtual labs

$$l_1 = 112.99$$

$$l_2 = 82.99$$

$$t_1 = 0.36 \text{ m/s}$$

$$t_2 = -1 \text{ m/s}$$

From virtual Labs

$$l_1 = 198.52$$

$$l_2 = 7.526$$

$$t_1 = 168.526$$

$$t_2 = -0.052$$

$$\therefore v = \frac{112.99 - 82.99}{0.36 - (-1)}$$

$$\therefore v = 198.52 - 7.526$$


$$168.52 - (-0.052)$$

$$v = 21.92 \text{ m/s}$$

$$v = 3.96 \text{ m/s}$$

## RESULTS :

We can understand that as we move the testing probe from the electrode, The velocity changes.  
Thus we can analyse the response of muscle.



## CONCLUSION :

Simulated the electromyogram Signal and analysed its characteristics.

## PRACTICAL NO: 4

Date : 06/02/2025.

TITLE : Analyzing image and EEG Signal frequency domain analysis using FFT in matlab.

AIM / OBJECTIVE: To apply FFT for EEG Signal analysis, compute the power spectrum, extract alpha, beta and theta waves, and visualize their temporal variations.

APPARATUS / TOOLS / EQUIPMENT / RESOURCES USED:

Apparatus

EEG Signal data.

Matlab.

CONCEPT / THEORY OF EXPERIMENT:

EEG Signal represents brain activity and contains frequency components linked to mental states. FFT converts time domain signal into freq domain enabling power spectrum computation and frequency band extraction.  
alpha 8-12 Hz, Beta 12-30 Hz, Theta 4-8 Hz.

PROCEDURE :

- 1 Load or generate EEG data in matlab.
- 2 Divide it into 1 sec frames.
- 3 Apply FFT to compute power spectrum.
- 4 Extract Alpha, Beta and Theta waves.
- 5 Plot their power variations over time.

```

% Generate a sample EEG signal
Fs = 1000; % Sampling frequency (Hz)
t = 0:1/Fs:10; % Time vector (10 seconds)
eeg_signal = sin(2*pi*10*t) + 0.5*sin(2*pi*20*t) + 0.2*sin(2*pi*5*t); % Sample EEG signal with frequencies 10 Hz, 20 Hz, and 5 Hz

% Frequency binning parameters
alpha_band = [8, 12]; % Alpha band frequencies (Hz)
beta_band = [13, 30]; % Beta band frequencies (Hz)
theta_band = [4, 7]; % Theta band frequencies (Hz)

% Prepare short 1-second frames of EEG signal
frame_length = Fs; % 1-second frame
num_frames = floor(length(eeg_signal) / frame_length);
frames = reshape(eeg_signal(1:num_frames*frame_length), frame_length, num_frames);

% Compute power spectrum and extract frequency bands for each frame
alpha_power = zeros(1, num_frames);
beta_power = zeros(1, num_frames);
theta_power = zeros(1, num_frames);

for i = 1:num_frames
    % Apply FFT to the frame
    frame_fft = fft(frames(:, i));

    % Compute power spectrum
    power_spectrum = abs(frame_fft).^2;


    % Extract frequency bins
    freqs = (0:frame_length-1) * Fs / frame_length;

    % Find indices corresponding to alpha, beta, and theta bands
    alpha_indices = find(freqs >= alpha_band(1) & freqs <= alpha_band(2));
    beta_indices = find(freqs >= beta_band(1) & freqs <= beta_band(2));
    theta_indices = find(freqs >= theta_band(1) & freqs <= theta_band(2));

    % Compute power in alpha, beta, and theta bands
    alpha_power(i) = sum(power_spectrum(alpha_indices));
    beta_power(i) = sum(power_spectrum(beta_indices));
    theta_power(i) = sum(power_spectrum(theta_indices));
end

% Plot the power in each frequency band over time
time = (1:num_frames) * (frame_length / Fs);
figure;
plot(time, alpha_power, 'r', 'LineWidth', 2);
hold on;
plot(time, beta_power, 'g', 'LineWidth', 2);
plot(time, theta_power, 'b', 'LineWidth', 2);
xlabel('Time (s)');
ylabel('Power');
title('Power in Alpha, Beta, and Theta Bands');
legend('Alpha', 'Beta', 'Theta');
grid on;

```



```

% Now let's calculate and plot the FFT of the Alpha, Beta, and Theta bands from the generated EEG signal
% Filter the EEG signal for alpha, beta, and theta bands
alpha_signal = zeros(size(eeg_signal));
beta_signal = zeros(size(eeg_signal));

```

```

theta_signal = zeros(size(eeg_signal));

for i = 1:num_frames
    % Extract the current frame
    frame = frames(:, i);

    % Apply bandpass filters to extract alpha, beta, and theta signals
    % Alpha band (8-12 Hz)
    alpha_signal((i-1)*frame_length + 1:i*frame_length) = bandpass(frame, alpha_band, Fs);

    % Beta band (13-30 Hz)
    beta_signal((i-1)*frame_length + 1:i*frame_length) = bandpass(frame, beta_band, Fs);

    % Theta band (4-7 Hz)
    theta_signal((i-1)*frame_length + 1:i*frame_length) = bandpass(frame, theta_band, Fs);
end

% Compute the FFT of the filtered signals
alpha_fft = fft(alpha_signal);
beta_fft = fft(beta_signal);
theta_fft = fft(theta_signal);


% Compute the frequency axis for the entire signal
N = length(eeg_signal); % Length of the full EEG signal
freqs = (0:N-1) * Fs / N; % Frequency axis

% Plot the FFT of alpha, beta, and theta signals
figure;
subplot(3, 1, 1);
plot(freqs, abs(alpha_fft));
title('FFT of Alpha Band (8-12 Hz)');
xlabel('Frequency (Hz)');
ylabel('Magnitude');
xlim([0 Fs/2]); % Limit x-axis to half the sampling frequency (Nyquist limit)

subplot(3, 1, 2);
plot(freqs, abs(beta_fft));
title('FFT of Beta Band (13-30 Hz)');
xlabel('Frequency (Hz)');
ylabel('Magnitude');
xlim([0 Fs/2]); % Limit x-axis to half the sampling frequency (Nyquist limit)

subplot(3, 1, 3);
plot(freqs, abs(theta_fft));
title('FFT of Theta Band (4-7 Hz)');
xlabel('Frequency (Hz)');
ylabel('Magnitude');
xlim([0 Fs/2]); % Limit x-axis to half the sampling frequency (Nyquist limit)

```



```

% If EEG csv file as input is given ,
% Read EEG signal from CSV file
eeg_data = csvread('C:\Users\student.SSPU\Downloads\Eeg_signal.csv'); % Assuming the CSV file contains one column representing

% Parameters
Fs = 1000; % Sampling frequency (Hz)
frame_length = Fs; % 1-second frame

% Frequency binning parameters
alpha_band = [8, 12]; % Alpha band frequencies (Hz)
beta_band = [13, 30]; % Beta band frequencies (Hz)
theta_band = [4, 7]; % Theta band frequencies (Hz)

```

```


'Prepare short 1-second frames of EEG signal
num_frames = floor(length(eeg_data) / frame_length);
frames = reshape(eeg_data(1:num_frames*frame_length), frame_length, num_frames);

; Compute power spectrum and extract frequency bands for each frame
alpha_power = zeros(1, num_frames);
beta_power = zeros(1, num_frames);
theta_power = zeros(1, num_frames);

for i = 1:num_frames
    % Apply FFT to the frame
    frame_fft = fft(frames(:, i));
    
    % Compute power spectrum
    power_spectrum = abs(frame_fft).^2;
    
    % Extract frequency bins
    freqs = (0:frame_length-1) * Fs / frame_length;
    
    % Find indices corresponding to alpha, beta, and theta bands
    alpha_indices = find(freqs >= alpha_band(1) & freqs <= alpha_band(2));
    beta_indices = find(freqs >= beta_band(1) & freqs <= beta_band(2));
    theta_indices = find(freqs >= theta_band(1) & freqs <= theta_band(2));
    
    % Compute power in alpha, beta, and theta bands
    alpha_power(i) = sum(power_spectrum(alpha_indices));
    beta_power(i) = sum(power_spectrum(beta_indices));
    theta_power(i) = sum(power_spectrum(theta_indices));
end

% Plot the power in each frequency band over time
time = (1:num_frames) * (frame_length / Fs);
figure;
plot(time, alpha_power, 'r', 'LineWidth', 2);
hold on;
plot(time, beta_power, 'g', 'LineWidth', 2);
plot(time, theta_power, 'b', 'LineWidth', 2);
xlabel('Time (s)');
ylabel('Power');
title('Power in Alpha, Beta, and Theta Bands');
legend('Alpha', 'Beta', 'Theta');
grid on;

```



```

% Read EEG signal from CSV file
eeg_data = readtable('C:\Users\student.SSPU\Downloads\yash.csv');

Warning: Variable names were modified to make them valid MATLAB identifiers. The original names are saved in the VariableDescriptions property.

% Assuming the CSV file contains columns representing the EEG signals




% Display variable names to check modifications
disp(eeg_data.Properties.VariableNames);

'EEG_AF3'    'EEG_T7'    'EEG_Pz'    'EEG_T8'    'EEG_AF4'

% Convert table to array by selecting the relevant columns
eeg_signals = eeg_data(:, 2:end); % Adjust indices as necessary

% Parameters
Fs = 1000; % Sampling frequency (Hz)
frame_length = Fs; % 1-second frame

```



C

## CALCULATIONS :

### \* Applications

- Neuroscience research
- Brain computer interfaces
- Mental health monitoring
- Sleep studies
- Cognitive Neuroscience
- Human factors Engineering
- Neurofeedback Therapy.

## RESULTS :

FFT successfully analyzed EEG signals revealing frequency characteristics and temporal variations.

This technique aids in neuroscience, psychology, and BCI research.

## CONCLUSION :

Hence, able to perform the equipment for analyzing EEG signals using FFT in matlab.

## PRACTICAL NO: 5

Date : 27/02/2025

TITLE: To Simulate Defibrillator

AIM / OBJECTIVE: 1. To simulate defibrillator output waveform  
2. To understand energy levels, necessary and application of defibrillator.

APPARATUS / TOOLS / EQUIPMENT / RESOURCES USED:

Virtual labs (Biomedical and signal processing)

Reference Notes.

CONCEPT / THEORY OF EXPERIMENT:

A defibrillator works by delivering a controlled electric shock to the heart restoring the normal rhythm during cardiac arrest. This shock interrupt chaotic electrical conductivity allowing the heart to restore its natural operating pattern. Defibrillators are crucial in saving lives by quickly restoring heart function in emergency situations.

PROCEDURE :

1. Run the defibrillator simulator and observe discharging waveform.
2. Run the defibrillator and observe energy delivered by changing voltage.
3. Note down the waveforms of the defibrillator by changing the values
4. Understand the importance of defibrillator.

### Defibrillator

College of Engineering  
(An Autonomous Institute of Government of Maharashtra)


Resistance(ohm):

Amplitude(v):

Types of Stimulus :

Energy : 134.02125 Joule

Energy =  $\frac{1}{2} CV^2$   
C = 33 microfarad



### Defibrillator

College of Engineering  
(An Autonomous Institute of Government of Maharashtra)


Resistance(ohm):

Amplitude(v):

Types of Stimulus :

Energy : 412.5 Joule

Energy =  $\frac{1}{2} CV^2$   
C = 33 microfarad



### Defibrillator

College of Engineering  
(An Autonomous Institute of Government of Maharashtra)


Resistance(ohm):

Amplitude(v):

Types of Stimulus :

Energy : 412.5 Joule


Energy =  $\frac{1}{2} CV^2$   
C = 33 microfarad



## OBSERVATIONS

- Energy levels :- Defibrillators typically generate energy levels ranging from 150 to 360 joules, crucial for restoring normal heart rhythm during cardiac arrest.
- Necessity / Application :- They are necessary in emergency medical settings to revive a stopped heart and sustain it by delivering controlled electric shock.
- Controls :- Common controls associated with defibrillators include energy level selection, charging mechanism and shock delivering buttons, enabling health care providers to administer treatment effectively.
- Configuration and types :- Defibrillators come in various configurations such as manual, semi-automatic and fully automatic, each offering different levels of user intervention, based on clinical requirement. Types include external, implanted and wearable devices, serving diverse patient needs and medical scenarios with precision and efficiency.
- Energy stored in capacitor follows:

$$E = \frac{Cu^2}{2}$$



## **CALCULATIONS :**

## RESULTS :

We successfully simulated the output waveform of defibrillator.

#### **CONCLUSION :**

From we were able to simulate defibrillator output waveform.

## PRACTICAL NO: 6

Date : 06/03/2025

TITLE: To simulate Pacemaker

AIM / OBJECTIVE: To simulate the pacemaker output make.

APPARATUS / TOOLS / EQUIPMENT / RESOURCES USED:

VLAB (Biomedical and Signal Processing)



CONCEPT / THEORY OF EXPERIMENT:

A pacemaker is a device implanted in the body to regulate abnormal heart rhythm. It consists of a pulse generator and leads that deliver electrical impulses to the heart to maintain a normal rhythm. This device continuously monitors the heart electrical activity and delivers pulses. When needed to keep the heart beating at right pace.

PROCEDURE :

- 1) Run the pacemaker simulator and observe the monophasic and biphasic stimulus pulse.
- 2) Run the pacemaker simulator and observe stimulus pulse by varying the pulse width and voltage.
- 3) Run pacemaker simulator and observe the energy delivered from stimulus pulse.

Name: Yash Karande | PRN: 2200601104 | Synchronous Pacemaker



## OBSERVATIONS

Pacemakers can be classified into two main types based on its location and how they are implanted.

### \* External pacemakers.

1. External pacemakers are temporary devices typically used in emergency situations.
2. They are not surgically implanted but rather attached adhesive electrode placed on the chest.
3. These are connected to the heart through temporary leads, which delivers electrical impulses to regulate the heart.
4. These are portable and provide short term support until a permanent solution is implemented.

### \* Internal Pacemakers.

1. Internal Pacemakers also known as implantable pacemakers are permanent devices implanted inside the body.
2. They consists of small generator and one or more leads that are inserted into the hearts through veins.
3. Internal pacemakers continuously monitor the heart rhythm and deliver electric pulse when necessary to regulate heartbeat.
4. Those pacemakers are programmable and can be adjusted by health care professional to meet the specific needs of the patients.

## CALCULATIONS :

Given values

Initial charge = 100%  
Final charge = 10%  
Rate of discharge = 10% per minute  
Time taken = ?

Using the formula:

$$t = \frac{100 - 10}{10} = 9$$

Time taken = 9 minutes

## RESULTS :

Hence we were able to simulate pacemaker output waveform.

## CONCLUSION :

Hence we were able to understand concept of pacemakers and simulate its output waveform.

## PRACTICAL NO: 7

Date : \_\_\_\_\_

TITLE: To simulate overall functionality of Hemodialysis Machine  
(Artificial Kidney)

AIM / OBJECTIVE: 1. To understand the block schematic/ module involved in Hemodialysis machine.  
2. To understand various measurement and control involved in Hemodialysis machine.

APPARATUS / TOOLS / EQUIPMENT / RESOURCES USED:

Virtual Lab.

CONCEPT / THEORY OF EXPERIMENT:

Hemo is greek word for blood. Dialysis means a filtering process. Hemodialysis means the process of filtering blood in the body. Hemodialysis is a method for removing waste products such as creatine and urea, as well as free water from the blood when the kidneys are in filtration failure.


PROCEDURE :

- 1) Set blood rate in normal given range.
- 2) Set control for ultrafiltration goal and ultrafiltration rate
- 3) Set dialysate flow rate in normal given range.
- 4) Run virtual hemodialysis machine
- 5) Simulate Blood leak detector.
- 6) Simulate Air Bubble detector.

# Virtual Haemodialysis Machine

**Virtual Haemodialysis Machine**

**College of Engineering, Pune**  
(An Autonomous Institute of Government of Maharashtra)



## OBSERVATIONS

### \* Blood flow.

The blood pump takes and return the blood from the patient via the arterial and venous needles respectively. The blood is confined to the disposable plastic tubing and does not come in contact with any part of the machine. The blood coming from the pump flow to the dialyzer returns to the patient through the venous needle. The pump speed, and the resulting blood flow rate is adjustable.

### \* Blood Pressure.

The blood pressure is measured at both when it taken out from the limb and also when it is returned.

### \* Blood Flow Rate.

The effect of blood flow is easy to follow any change in pump speed is immediately reflected the pressure displays. The higher the flow the higher the pressure.

### \* Dialysis Process.

1. Concerning the blood, dialysis performs 2 different function that are normally done by healthy kidney.
2. Removing excess fluid.
3. Removing waste like urea and excess electrolytes like potassium, magnesium, sodium etc.

## CALCULATIONS :

and the total volume of the dialysis bath is  $100 \text{ ml}$ .  
The initial concentration of urea in the blood is  $100 \text{ mg/dl}$ .  
The initial concentration of urea in the dialysis bath is  $10 \text{ mg/dl}$ .  
The flow rate of blood is  $200 \text{ ml/min}$ .  
The flow rate of dialysis bath is  $100 \text{ ml/min}$ .  
The time taken for the urea to be removed from the blood is  $10 \text{ min}$ .  
The final concentration of urea in the blood is  $50 \text{ mg/dl}$ .  
The final concentration of urea in the dialysis bath is  $50 \text{ mg/dl}$ .

## RESULTS :

Hence we are able to simulate the hemodialysis machine.  
(Artificial kidney)

## CONCLUSION :

Hence we are able to understand the concept of hemodialysis machine.